Author's personal copy

Journal of Chromatography A, 1217 (2010) 2674-2684



Contents lists available at ScienceDirect

Journal of Chromatography A

journal homepage: www.elsevier.com/locate/chroma



Review

Sorbent-based sampling methods for volatile and semi-volatile organic compounds in air

Part 1: Sorbent-based air monitoring options

Elizabeth Woolfenden*

Markes International Ltd., Gwaun Elai Campus, Llantrisant RCT, CF72 8XL, UK

ARTICLE INFO

Article history: Available online 11 January 201

Available online 11 January 2010

Keywords:

Air monitoring
Sorbent tubes
Sorbent traps
VOCs
SVOCs
Volatile organic chemicals
Semi-volatile organic chemicals
Thermal desorption

ABSTRACT

Sorbent tubes/traps are widely used in combination with gas chromatographic (GC) analytical methods to monitor the vapour-phase fraction of organic compounds in air. Target compounds range in volatility from acetylene and freons to phthalates and PCBs and include apolar, polar and reactive species. Airborne vapour concentrations will vary depending on the nature of the location, nearby pollution sources, weather conditions, etc. Levels can range from low percent concentrations in stack and vent emissions to low part per trillion (ppt) levels in ultra-clean outdoor locations. Hundreds, even thousands of different compounds may be present in any given atmosphere. GC is commonly used in combination with mass spectrometry (MS) detection especially for environmental monitoring or for screening uncharacterised workplace atmospheres. Given the complexity and variability of organic vapours in air, no one sampling approach suits every monitoring scenario. A variety of different sampling strategies and sorbent media have been developed to address specific applications. Key sorbent-based examples include: active (pumped) sampling onto tubes packed with one or more sorbents held at ambient temperature; diffusive (passive) sampling onto sorbent tubes/cartridges; on-line sampling of air/gas streams into cooled sorbent traps; and transfer of air samples from containers (canisters, Tedlar® bags, etc.) into cooled sorbent focusing traps. Whichever sampling approach is selected, subsequent analysis almost always involves either solvent extraction or thermal desorption (TD) prior to GC(/MS) analysis. The overall performance of the air monitoring method will depend heavily on appropriate selection of key sampling and analytical parameters. This comprehensive review of air monitoring using sorbent tubes/traps is divided into 2 parts. (1) Sorbent-based air sampling option. (2) Sorbent selection and other aspects of optimizing sorbent-based air monitoring methods. The paper presents current state-of-the-art and recent developments in relevant areas such as sorbent research, sampler design, enhanced approaches to analytical quality assurance and on-tube derivatisation.

© 2009 Elsevier B.V. All rights reserved.

Contents

1.	Introduction	2675
2.	Summary of sorbent-based air sampling options	2676
	Pumped (active) sampling onto sorbent tubes	
	Diffusive (passive) sampling	
	4.1. Axial diffusive samplers	2678
	4.2. Radial diffusive samplers	2679
5.	Whole-air sampling into sorbent focusing traps	
	5.1. Online operation	2680
	5.2. Offline air sampling using containers (canisters, Tedlar bags, etc.)	2681
6.	A quick review of some alternative air sampling approaches.	
	6.1. Solid phase (micro-)extraction (SP(M)E)	2683
	6.2. Equilibrium/static headspace	

^{*} Tel.: +44 1443 230935; fax: +44 1443 231531.

E-mail address: ewoolfenden@markes.com.